



Weighing scoop

Designed for accurate transfer

Made from 3.3 Low Expansion Borosilicate Glass



Borosil Scientific's weighing scoops are designed to offer precise, contamination-free transfer of powders and granules directly to your balance. With smooth surfaces for easy cleaning and ergonomic shapes for comfortable use, they are an essential tool for any laboratory or production setup.

What are Weighing Scoop used for?

Borosil Scientific's weighing scoops are primarily used when small quantities of powder, liquid or semi-solid material need to be weighed. The powder can be introduced into the received vessel through the tabular stem.

Borosil Scientific's weighing scoops are engineered for reliable performance & are durable, easy to clean, resistant to most lab chemicals - perfect for routine and analytical lab work.

Primary Uses of Weighing Scoop:

- Weighing and Transferring Solids/ Liquids/ Semi Solids
- Minimizing Sample Loss
- Efficient Sample Handling
- Avoiding Cross-Contamination
- Measurement Support

Weighing Scoop have open tabular arm.
Ideal for weighing small quantities drugs and powders



Product Code	Capacity ml	Quantity Per case
7200P45	0.5	10
7200001	1	10
7200002	2	10
7200003	3	10
7200076	6	10
7200006	10	10
7205706	10	10

BOROSIL®
Scientific
Borosil Scientific Limited

CIN : U74999MH1991PLC061851

Visit Website



Follow us on LinkedIn



Registered & Corporate Office :

1101, Crescenzo, G-Block, Opp. MCA Club, Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India.

T +91 22 6740 6300 ● F +91 22 6740 6514

E bsl@borosil.com ● W www.borosilscientific.com

Write to us for export enquiries at bsl.exports@borosil.com

MUMBAI

mumbaisales@borosil.com

+91 22 6740 6400
+91 22 6740 6444

DELHI

delhi@borosil.com

+91 11 4150 5893

CHENNAI

chennai@borosil.com

+91 44 2822 6012 / 13
+91 44 2822 6014

KOLKATA

calcutta@borosil.com

+91 33 2229 9166

All rights reserved. Borosil Scientific Limited. 2025.

Any unauthorized use, copying, or distribution of this brochure is strictly prohibited and constitutes a legal violation. Borosil Scientific Limited reserves all rights to take legal action against individuals or entities who violate these legal rights.