




Microscope slides & Cover glasses

Reliable Clarity for Every Specimen

Borosil Scientific's microscope slides are engineered for clarity and consistency, making them ideal for precise sample examination. Crafted from high-quality glass, they ensure smooth edges, excellent transparency, and compatibility with all standard microscopes.

Uses of Microscope Slides


Microscope slides are essential tools for laboratory use in preparing and examining specimens under a microscope. They are widely used in biological, chemical, and medical laboratories to hold tissue samples, cells, microorganisms, or chemical substances for observation & analysis.



Material

Made of Sodalime glass with thickness of 1.0 - 1.2mm

Types



Packaging

Packed in a paper box containing silica desiccant to absorb moisture 50 slides per packet & 25 packets per case

Product Code	Description	Dimensions (mm) 18 X 18	Quantity Per Packet (No)	Packets / Case
* 9100P02	Plain, Ground Edges	76 X 26 X 1	50	25
9100D01	Double Frosted Ground Edges	76 X 26 X 1	50	25


* ISI marking is applicable only for plain packaging with grounded edges.



Scan for Product Details

Why Choose Borosil Scientific Cover Slips?

- High Optical Clarity - Ideal for detailed microscopic observations.
- Uniform Thickness - Ensures consistent focus and imaging.
- Cleanroom Packing - Reduces contamination risks.
- Shatter Resistance – Carefully annealed to reduce breakage during handling.



Square - 9115501



Square - 9115S02



Square - 9115S03



Rectangular - 9115R01



Rectangular - 9115R02



Rectangular - 9115R03



Conformance to Standard ASTM E211

Product Code	Description	Dimensions (mm)	Quantity / Packet	Packets / Case
9115S01	Square	18 X 18	100	50
9115S02		22 X 22		
9115S03		24 X 24		
9115R01	Rectangular	22 X 40	100	50
9115R02		22 X 50		
9115R03		22 X 60		



Scan for Product Details



CIN : U74999MH1991PLC061851

Visit
Website



Follow us on
LinkedIn



Registered & Corporate Office :

1101, Crescenzo, G-Block, Opp. MCA Club, Bandra Kurla Complex,
Bandra (E), Mumbai - 400 051, India.

T +91 22 6740 6300 ● F +91 22 6740 6514

E bsl@borosil.com ● W www.borosilscientific.com

Write to us for export enquiries at bsl.exports@borosil.com

MUMBAI

mumbaisales@borosil.com

+91 22 6740 6400

+91 22 6740 6444

DELHI

delhi@borosil.com

+91 11 4150 5893

CHENNAI

chennai@borosil.com

+91 44 2822 6012 / 13

+91 44 2822 6014

KOLKATA

calcutta@borosil.com

+91 33 2229 9166

All rights reserved. Borosil Scientific Limited. 2025.

Any unauthorized use, copying, or distribution of this brochure is strictly prohibited and constitutes a legal violation. Borosil Scientific Limited reserves all rights to take legal action against individuals or entities who violate these legal rights.